

HIV Rapid

Synonyms	Rapid HIV
Abbreviation	HIV 1/2
Specimen type	Preferred: Serum Acceptable: Plasma
Container	Red, Red/Yellow (SST) Green (PST), Lavender, Blue
Volume	1 mL
Collection requirements	Deliver to Lab STAT Confirmatory testing done if Pos.
S/H instructions	Refrigerate
Stability	Centrifuged and separated: 3 days refrigerated Freeze if stored longer
Reference range	Non-Reactive
Routine / STAT T.A.T	Within 45 min / 45 min
Section	Serology
LAB code	08870