

## Strep Pneumo Ag

<b>Synonyms</b>	Pneumococcal Ag Urine
<b>Abbreviation</b>	Strep pneumo
<b>Specimen type</b>	Urine <b>Acceptable:</b> CSF
<b>Container</b>	Sterile, leakproof
<b>Volume</b>	1 mL Minimum: 1 mL
<b>Collection requirements</b>	Collect voided urine or urine via catheter. Standard CSF collection into CSF tube.
<b>S/H instructions</b>	Deliver to Lab within 30 minutes; assay within 24 hours.
<b>Stability</b>	Urine: up to 14 days refrigerated or frozen. CSF: up to 7 days refrigerated or frozen.
<b>Reference range</b>	No Strep Pneumo antigen detected
<b>Routine / STAT T.A.T</b>	Within 4 hours / 1 hour
<b>Section</b>	Serology
<b>LAB code</b>	07119