

Culture Anaerobic

Synonyms	ANA/AER Culture, Anaerobic culture
Abbreviation	ANA/AER CX
Specimen type	Fluid aspirate/tissue
Container	Sterile syringe/container Anaerobic swab
Volume	≤ 0.1 mL
Collection requirements	Aspirate or tissue preferred. Anaerobic swabs acceptable and available in Lab.
S/H instructions	Deliver to Lab within 30 min. Includes Routine culture and Gram stain.
Stability	
Reference range	
Routine / STAT T.A.T	Prelim within 24 hrs
Section	Micro
LAB code	04252