

Digoxin Level

Synonyms	Dig
Abbreviation	Dig
Specimen type	Plasma Serum is acceptable
Container	Preferred: Green (PST) Acceptable: Red, Red/Yellow (SST) Lavender *Tube types should NOT be used interchangeably.
Volume	1 mL Minimum: 0.5 mL
Collection requirements	Draw levels 6-8 hours after daily dose or immediately before next dose. Verify last dose date/time.
S/H instructions	Refrig
Stability	Centrifuged: 48 hrs at room temp or refrigerated 6 months frozen
Reference range	0.8 – 2.0 ng/mL (Therapeutic) ≥ 3.0 ng/mL (Panic; when ≥6hrs after dose)
Routine / STAT T.A.T	Within 4 hours / 1 hour
Section	Chem
LAB code	01050