

Eosinophil Urine Smear

Synonyms	Urine Eos Smear Eos Urine
Abbreviation	EOS Urine
Specimen type	Preferred: Fresh Urine, < 4 hours old
Container	Sterile cup
Volume	10 mL Minimum: 10 mL
Collection requirements	Sample should be less than 4 hour old.
S/H instructions	
Stability	4 hours
Reference range	No Eosinophils present
Routine / STAT T.A.T	Within 4 hours / 1 hour
Section	Heme
LAB code	LAB02106