

FFP Units available

Synonyms	
Abbreviation	FFP units
Specimen type	Whole blood
Container	Pink top tube Acceptable: Lavender, Red
Volume	6 mL Minimum: 3 mL
Collection requirements	Avoid hemolysis. ABO compatible units should be given. ID patient and use Hollister ID band. *Specimen needed if patient not previously typed.
Special instructions	Thaw plasma in ziplock bag in water bath @ 30-37 C, with gentle aggitation every 5 min. After thawing, store in BB frig @ 1-6 C. Unit must be given within 24hrs.
Stability	Keep FFP refrigerated after thawing; good for up to 24 hrs.
Reference range	NA
Routine / STAT T.A.T	Within 2 hours / 1 hour

Section	Blood bank
LAB code	03092