

Mononucleosis Screen

Synonyms	Mono Screen, Monospot, Mono test, Heterophile Ab, IM antibody, Infectious Mono Ab
Abbreviation	Mono
Specimen type	Preferred: EDTA Plasma Acceptable: Serum
Container	Lavender (EDTA), Red, Red/Yellow (SST)
Volume	0.5 mL Minimum: 0.2 mL
Collection requirements	Bring sample directly to Lab. Must be centrifuged and serum separated within 2 hours of collection.
S/H instructions	Refrigerate
Stability	Centrifuged and separated: 72 hours
Reference range	Negative
Routine / STAT T.A.T	Within 4 hours / 1 hour
Section	Serology
LAB code	07040

