

Rapid Bacterial Vaginosis

| | |
|--------------------------------|--|
| Synonyms | Rapid B.V., B.V. Rapid, Bacterial Vaginosis Rapid |
| Abbreviation | Rapid B.V. |
| Specimen type | Swab of lower third of vaginal wall. |
| Container | BV Blue swab from OSOM collection kit |
| Volume | |
| Collection requirements | Collect from the lower third of the vaginal wall using the swab from the kit; <i>avoid cervical collections.</i> Avoid vaginal creams, lubricants, douches, etc. within the last 72 hrs. |
| S/H instructions | Deliver to Lab immediately. |
| Stability | Process ASAP; may be stored up to 48 hr refrigerated. |
| Reference range | Negative |
| Routine / STAT T.A.T | Within 4 hours / 1 hour |
| Section | Serology |

| | |
|-----------------|-------|
| LAB code | 90304 |
|-----------------|-------|