

Rapid Trichomonas

Synonyms	Trichomonas, Trich rapid
Abbreviation	Rapid Trich
Specimen type	Swab of vaginal wall
Container	Sterile rayon swab in collection kit. Cotton-tip swabs are <i>not</i> recommended.
Volume	
Collection requirements	Collect from vaginal cavity with sterile rayon swab provided in kit; avoid vaginal lubricants.
S/H instructions	Deliver to Lab asap.
Stability	Up to 24 hours at room temp.
Reference range	Negative
Routine / STAT T.A.T	Within 4 hours / 1 hour
Section	Serology
LAB code	04292