

Triglycerides

Synonyms	Trig
Abbreviation	TRIG
Specimen type	Preferred: Plasma Acceptable: Serum
Container	Green/Yellow (pst), Dark Green, Green/White, Green Microtainer; Red/Yellow (SST) or Plain Red Top or Red Top Microtainer.
Volume	1 mL Minimum: 0.5 mL Plasma/Serum
Collection requirements	Fasting specimen preferred.
S/H instructions	
Stability	Centrifuged: 2 days at room temperature. 7 days refrigerated. 1 year frozen.
Reference range	<150 mg/dL Desirable 150-199 mg/dL Borderline High 200-499 mg/dL High >500 mg/dL Very High
Routine / STAT T.A.T	Within 4 hours / 1 hour

Section	CHEM
LAB code	01690