

### Type + Screen Panel

<b>Synonyms</b>	ABO/RH + Screen, T and S
<b>Abbreviation</b>	TANDS
<b>Specimen type</b>	Whole blood <b>Preferred:</b> Plasma & cells from Pink top <b>Acceptable:</b> Plasma & cells from Lavender Serum & cells from Red top
<b>Container</b>	1 Pink + 1 Red Acceptable: Lavender
<b>Volume</b>	6 mL Minimum: 3 mL
<b>Collection requirements</b>	Avoid hemolysis
<b>S/H instructions</b>	Refrigerate
<b>Stability</b>	7 days at 1-8 C
<b>Reference range</b>	NA
<b>Routine / STAT T.A.T</b>	Within 4 hours / 1 hour
<b>Section</b>	Blood Bank
<b>LAB code</b>	03029