

## Legionella Ag Urine

|                                |                                                                                                                |
|--------------------------------|----------------------------------------------------------------------------------------------------------------|
| <b>Synonyms</b>                | Urine Legionella Ag, Legionella urine Ag                                                                       |
| <b>Abbreviation</b>            | Legionella                                                                                                     |
| <b>Specimen type</b>           | 5.0 ml Random Urine                                                                                            |
| <b>Container</b>               | Sterile Urine Container                                                                                        |
| <b>Volume</b>                  | Minimum: 0.5 ml                                                                                                |
| <b>Collection requirements</b> | Collect first catch urine in sterile container.<br>Avoid diuretic urine.<br>Deliver immediately to laboratory. |
| <b>S/H instructions</b>        | Hold and transport refrigerated.                                                                               |
| <b>Stability</b>               | 24 hours at room temp.<br>Up to 14 days refrigerated or frozen.                                                |
| <b>Reference range</b>         | Negative                                                                                                       |
| <b>Routine / STAT T.A.T</b>    | Within 4 hours / 1 hour                                                                                        |
| <b>Section</b>                 | Serology                                                                                                       |
| <b>LAB code</b>                | 06351                                                                                                          |