

Urinalysis Routine

Synonyms	UA, URINALYSIS WITH REFLEX MICRO
Abbreviation	UA
Specimen type	Urine
Container	Sterile Specimen Container
Volume	10 mL Minimum: 3 mL Urine
Collection requirements	Send to Lab within 2 hours of collection. Microscopic will only be done if any of the following are positive: Glucose, Protein(>trace), Blood, Nitrite, Leukocyte
S/H instructions	Refrigerate
Stability	2 hours at room temp. 12 hours Refrigerated.
Reference range	Color: Yellow Glu Neg Ket Neg Bili Neg pH 5.0-8.0 Bld Neg Sp Gr 1.003-1.030 Pro Neg Urobilinogen <2.0 mg/dL WBC Neg Nit Neg

Routine / STAT T.A.T	Within 4 hours / 1 hour
Section	UA
LAB code	07530