

Potassium Urine Timed

Synonyms	K urine timed, Timed Urine potassium, Urine potassium timed
Abbreviation	URK 24hr
Specimen type	24 Hour Urine
Container	24 Hour Urine Container Without Preservatives
Volume	10mL aliquot or 24hr collection Minimum: 0.5 mL
Collection requirements	24 Hour Urine container available in Lab. Refrigerate or place on ice during collection.
S/H instructions	Refrigerate
Stability	45 days at room temperature 2 months refrigerated 1 year frozen
Reference range	25-125 mmol/24 hours
Routine / STAT T.A.T	Within 4 hours / 1 hour
Section	CHEM
LAB code	01620

