

Alkaline Phosphatase

Synonyms	Alk phos, Alkp																					
Abbreviation	Alk phos																					
Specimen type	Plasma																					
Container	Green (PST)																					
Volume	1.0 mL Minimum: 0.5 mL																					
Collection requirements	Avoid hemolysis. Remove from cells within 4 hours of collection.																					
S/H instructions	Refrigerate																					
Stability	Centrifuged and removed from cells: 48 hours refrigerated																					
Reference range	<table border="1"><thead><tr><th></th><th>Male</th><th>Female</th></tr><tr><th>Age</th><th>Range U/L</th><th>Range U/L</th></tr></thead><tbody><tr><td>0 – 11yr</td><td>135 - 530</td><td>130 – 560</td></tr><tr><td>11 – 13yr</td><td>200 – 495</td><td>105 – 420</td></tr><tr><td>13 – 15yr</td><td>130 – 525</td><td>70 -230</td></tr><tr><td>15 – 19yr</td><td>65 – 260</td><td>50 -130</td></tr><tr><td>>19yr</td><td>38 - 126</td><td>38 - 126</td></tr></tbody></table>		Male	Female	Age	Range U/L	Range U/L	0 – 11yr	135 - 530	130 – 560	11 – 13yr	200 – 495	105 – 420	13 – 15yr	130 – 525	70 -230	15 – 19yr	65 – 260	50 -130	>19yr	38 - 126	38 - 126
	Male	Female																				
Age	Range U/L	Range U/L																				
0 – 11yr	135 - 530	130 – 560																				
11 – 13yr	200 – 495	105 – 420																				
13 – 15yr	130 – 525	70 -230																				
15 – 19yr	65 – 260	50 -130																				
>19yr	38 - 126	38 - 126																				
Routine / STAT T.A.T	Within 4 hours / 1 hour																					
Section	Chem																					
LAB code	01220																					