

Crossmatch RBC

Synonyms	XMATCH
Abbreviation	XMAT
Specimen type	Whole blood Preferred: Plasma & cells from Pink top Acceptable: Plasma & cells from Lavender Serum & cells from Red top
Container	1 Pink Acceptable: Lavender, Red
Volume	6 mL Minimum: 3 mL
Collection requirements	Avoid hemolysis
S/H instructions	Refrigerate
Stability	7 days at 1-8 C
Reference range	NA
Routine / STAT T.A.T	Within 4 hours / 1 hour
Section	Blood bank
LAB code	03122